

Alt



CLEANING BRUSHES

ECOFLUSHCAPS

JET CHANNEL ADAPTER

ACCESSORIES

CLEANING BRUSHES

Safe and gentle cleaning

- Rounded brush tips
- Resistant to bending
- ✓ Latex-free



For flexible Endoscopes



Cleaning brushes

Safe and gentle cleaning

Features

- igsirial Ball heads on the brush ends and rounded brush tips for gentle cleaning
- \bigtriangledown Brushes made of soft, durable nylon
- \bigtriangledown Ideal brush head size, with optimal density
- igsilon Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- \bigtriangledown Bend resistant thanks for first-class catheter rigidity

SINGLE/DUAL-ENDED CLEANING BRUSHES

							Features	
							➡ Efficient clea brushes on c	ning thanks to two ne shaft
Colour	Cath	neter	Bru	ısh 1	Bru	sh 2	Pieces per box	Manufacturer Product
	Ø	Length	Ø	Length	Ø	Length		Code
		FOR EN	DOSCOPES		NG CHANNEL	S FROM 2.8 -	- 4.2 MM	
Yellow	1.7 mm	230 cm	6 mm	20 mm	-	-	100	SCBS-23-20-06-230
Blue	1.7 mm	230 cm	6 mm	20 mm	6 mm	20 mm	100	DCBS-23-20-06-230
Yellow	1.8 mm	180 cm	5 mm	20 mm	5 mm	20 mm	100	DCBS-23-20-05-180
Yellow	1.8 mm	230 cm	5 mm	20 mm	5 mm	20 mm	100	DCBS-23-20-05-230
		FOR PEDIATR	IC ENDOSCO	OPES WITH W	ORKING CHA	NNELS FROM	1 2.0 – 2.8 MM	
Green	1.5 mm	120 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-120
Green	1.5 mm	180 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-180
Green	1.5 mm	230 cm	3 mm	20 mm	3 mm	20 mm	100	DCBS-18-20-03-230
	FC	OR PEDIATRIC	BRONCHOS	COPES WITH	WORKING C	HANNELS FR	OM 1.2 – 1.8 MM	
Yellow	1.0 mm	120 cm	2 mm	15 mm	2 mm	15 mm	100	DCBS-15-15-02-120
		FOR ENT	EROSCOPES	WITH WORK	ING CHANNE	LS FROM 2.8	– 4.2 MM	
Blue	1.7 mm	300 cm	6 mm	20 mm	6 mm	20 mm	100	DCBS-23-20-06-300
		FOR EN	DOSCOPES	WITH WORKI	NG CHANNEL	S FROM 1.0 -	- 1.5 MM	
Blue	0.8 mm	120 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	DCBS-08-10-15-120
Blue	0.8 mm	180 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	DCBS-08-10-15-180
Blue	0.8 mm	230 cm	1.5 mm	10 mm	1.5 mm	10 mm	100	DCBS-08-10-15-230



For flexible endoscopes




Cleaning brushes

Safe and gentle cleaning



- \bigtriangledown Ball heads on the brush ends and rounded brush tips for gentle cleaning
- \bigtriangledown Brushes made of soft, durable nylon
- igsquare Ideal brush head size and density
- igsilon Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- \bigtriangledown Bend resistant thanks for first-class catheter rigidity

MULTIPOINT DUAL-ENDED CLEANING BRUSH





	 -	r-
9	 -	



COMBINATION CLEANING BRUSHES

Flexible and efficient cleaning of working channels and valves

- High cleaning efficiency
- Rounded brush tips
- ✓ Latex-free



For flexible endoscopes



mdd


Combination cleaning brushes

Flexible and efficient cleaning of working channels and valves


Features

- \bigtriangledown Ball heads on the brush ends and rounded brush tips for gentle cleaning
- \bigtriangledown Brushes made of soft, durable nylon
- igsquare Ideal brush head size and density
- igsilon Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- \bigtriangledown Bend resistant thanks for first-class catheter rigidity

WITH HANDLE - TWO-SIDED



WITHOUT HANDLE – TWO-SIDED





For flexible endoscopes




Combination cleaning brushes

Flexible and efficient cleaning of working channels and valves


Features

- \bigtriangledown Ball heads on the brush ends and rounded brush tips for gentle cleaning
- \bigtriangledown Brushes made of soft, durable nylon
- igsquare Ideal brush head size and density
- igsilon Brush-head diameter allows for maximum mechanical brushing of the endoscope channel
- \bigtriangledown Bend resistant thanks for first-class catheter rigidity

WITH SWEEPER



TWO-SIDED (CHRISTMAS-TREE SHAPE)






mdd

CLEANING BRUSHES

For flexible endoscopes



Combination cleaning brushes Flexible and efficient cleaning of working channels and valves TWO-SIDED (CHRISTMAS-TREE SHAPE) WITH SWEEPER



VALVE/CHANNEL CLEANING BRUSHES

Handy cleaning of working channels and valves

- High cleaning efficiency
- Rounded brush tips
- Maximum brush density



b

For flexible endoscopes




Valve/channel brushes

Handy cleaning of working channels and valves

Features

- \bigtriangledown Rounded brush tips for gentle cleaning
- \bigtriangledown Brushes made of soft, durable nylon
- \bigtriangledown Rounded handle

TWO-SIDED VALVE/CHANNEL CLEANING BRUSH



VALVE CLEANING BRUSH WITH HANDLE



CLEANING BRUSH FOR ALBARRÀN LEVER ON DUODENOSCOPES

Colour	Length	Channe	el brush	Pieces per box	Manufacturer Product	
		Ø	Length		Code	
Red	8 cm	2 mm	15 mm	200	SES-30-2-15	



Practical sets for efficient cleaning of working channels and valves


- High cleaning efficiency
- Rounded brush tips
- Ergonomic handle



CLEANING BRUSH SETS

For flexible endoscopes





MULTIPOINT DUAL-ENDED CLEANING BRUSH					DUAL	-ENDED VALVE	E CLEANING B	RUSH	Pieces per box	Manufacturer Product Code
Colour	Catheter		Brush 1 / 2 / 3 Valve brush 1		orush 1	Channel brush 2				
	Ø	Length	Ø	Length	Ø	Length	Ø	Length		
			FOR END	OSCOPES	WITH 2.8 -	4.2 MM W	ORKING CI	HANNELS		
Blue	1.7 mm	230 cm	6 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	MCBS-SET- 23-20-06-230
Blue	1.7 mm	260 cm	6 mm	20 mm	11 mm	35 mm	5 mm	20 mm	100	MCBS-SET- 23-20-06-260



Preparation of endoscopes with rinsing fluids

- Can be used for 24 hours
- ✓ With and without CO2 connection
- Olympus / Fujifilm / Pentax



(A

STERILE EO

With CO2 connection

CE

EcoFlushCap HybriX single use

Higher patient protection combined with lower costs

The innovative and single use **EcoFlushCap HybriX** supports the high demands on patient safety and significantly reduces the risk of cross-contamination. Connected to only one sterile water bottle system, the hybrid hose ensures safe flushing through the working channel and gentle flushing of the optics in the endoscope. The sterile, single use system is suitable for connection to all common flushing pump units through intelligent twin tubing. **EcoFlushCap HybriX** also supports CO2 insufflation thanks to the optional single use connection to CO2 insufflation devices.



Pieces per box		Manufacturer Product Code	
	For Olympus	For Fujifilm Series 500 / 600	For Pentax
15	WB-HYB-OL4	WB-HYB-FJ4	WB-HYB-PX6



With CO2 connection



EcoFlushCap single-use

Preparation of endoscopes with rinsing fluids

Features

- \bigtriangledown Can be used for 24 hours before cross contamination
- \bigtriangledown With CO2 hose system
- \bigtriangledown With hose set and screw cap for all common sterile water bottles



CLEAN CONNECTING TUBE FOR STERILE WATER BOTTLES WITH CO2 CONNECTION



CLEAN CONNECTING TUBE FOR STERILE WATER BOTTLES WITH 1.1 M CO2 CONNECTING HOSE





Without CO2 connection

 $\mathbf{C} \in (2)$

STERILE EO



EcoFlushCap single-use

Process-safe and efficient cleaning of lenses

Features

- \bigtriangledown Can be used for 24 hours before cross contamination
- \bigtriangledown With hose set and screw cap for all common sterile water bottles



CLEAN CONNECTING TUBE FOR STERILE WATER BOTTLES WITHOUT CO2 CONNECTION

Features

fluids

↓ With integrated non-return valve
 ↓ Individually sterile packed






Jet channel adapter & adapter





For Olympus







 \bigtriangledown Used together with the rinsing hose in connection with a rinsing pump for flushing with rinsing

For rinsing systems

m

dd	EcoFlushCap accessories
uu	For single-use clean connecting tube

C E ②

CO2 CONNECTION TUBE Features Q Length 120 cm Q 24 h use Q Single-use in order to support patient safety Luer Lock "female" Length Luer Lock "male" Interpretation Luer Lock "male" 10 Co2-TUI-210

IRRIGATION TUBE FOR OLYMPUS IRRIGATION PUMPS

Features		
 □ Can be used for 24 hours □ Connection for one-time use in the interest of patient safet □ For all common sterile water bottles □ Individually sterile packed 	y	
Product description	Pieces per box	Manufacturer Product Code
For Olympus rinsing pumps	15	IT-OUS-210

IRRIGATION TUBE WITH 2 LUER LOCK MALE AND CHECK VALVE

	Features		
Luer Lock "male"	 □ Can be used for 24 hours □ Connection for one-time use in the interest of patient safet □ For all common sterile water bottles □ Individually sterile packed 	/	
	Product description	Pieces per box	Manufacturer Product Code
	For rinsing pumps	15	IT-ETS-210



For rinsing systems



mdd

Bottle holder for sterile water bottles

Features For 1000 ml Ampuwa bottles and 500 ml Braun rinsing solution

- ♀ 2 hooks for rail mounting
- □ Internal diameter 90 mm

For 1000 ml Braun bottles

- ♀ 2 hooks for rail mounting
- □ Internal diameter 94 mm
- ➡ Height 170 mm

Product description	Pieces per box	Manufacturer Product Code
For 1000 ml Ampuwa bottles and 500 ml Braun rinsing solution	1	N1010364
For 1000 ml Braun bottles	1	N1010393

PLASTIC

8.8	Features		
	 For 1000 ml Ampuwa bottles and 500 ml Braun rinsing ♀ 1 hook for rail mounting ♀ With bracket for rear attachment ♀ 115 x 100 x 90 mm (L x W X H) 	g solution	
	Product description	Pieces per box	Manufacturer Product Code
	For 1000 and 500 ml Braun bottles	1	??







UpViser Aps Adelgade 45, 1 8660 Skanderborg, Denmark Tel.: +45 6171 7625 E-Mail: salg@upviser.com www.upviser.dk



Exclusive Partner Europe Body Products GmbH

mdd Company GmbH Rengglochstr. 38 CH-6012 Obernau, Switzerland

Tel: +39 348 0026970 E-Mail: gb@mddcompany-gmbh.com www.mddcompany-gmbh.com